Smithsonite - ZnCO₃

Hardness: 4 - 4.5. Cleavage: Very good. Insoluble in water. Rapidly soluble in acids except in cold saturated citric acid in which it is only slowly soluble. Avoid strong detergents. Soluble in alkalis but insoluble in NH₃. Derust with any dithionite recipe followed by a soak in NH₃. In one sample treatment with Jacquard Hydro formed a black coating which faded after sonication and thorough drying. Stable to acetone or 99% rubbing alcohol.

References

Mindat: https://www.mindat.org/min-3688.html

Handbook of Mineralogy: https://www.handbookofmineralogy.org/pdfs/smithsonite.pdf

Hardinger, S. (2024) Selective removal of carbonates: A preliminary study.

www.mindat.org/a/hardinger005

Hardinger, S, (2025) Mineral Specimen Cleaning and Development for the Amateur, 339 p.

Rohner, T. (2000) Properly clean minerals online cleaning book.

www.strahlen.org/stepbystep/mineralien-reinigung2.pdf. In German.

Sinkankas, J. (1972) Gemstone & Mineral Data Book, 346 p. Winchester Press, New York.

Voelter, J. (2010) Cleaning smithsonite, in Smithsonite Think Zinc!, ExtraLAPIS 13, 100 p. Lithographie, LLC, East Hampton, Connecticut.

Weast, R., Ed. (1982) CRC Handbook of Chemistry and Physics, 2380 p., CRC Press, Inc., Cleveland.